

## BRACE APPLICATION - INSTRUCTIONS

### 1. Initial preparation

Sit on the edge of a stable chair and flex the knee at a 45° angle (**Figure A**).



Figure A

Tie the straps to themselves so they do not bother you during installation (**Figure B**).



Figure B

### 2. Alignment of the EVOKE™

Align the center of the hinges (**Figure C**) at the center of the condyles of your knee (**Figure D**).



Figure C



Figure D

### 3. Strap #1

Attach the strap #1 proximal posterior tibial over the calf. It will keep the brace in place (**Figure E**).



Figure E

### 4. Strap #2 and #3

Attach the strap #2 anterior tibia, if present (**Figure F**).



Figure F



Figure G

Attach the strap #3 posterior distal tibia, if present (**Figure G**).

### 5. Strap #4, #5 and #6

Attach the strap #4 posterior proximal femoral and strap #5 posterior distal femoral (**Figure H**).




Figure H



Figure I

Attach the strap #6 anterior femoral, if present (**Figure I**).

### 6. Last adjustment

In order to check the fit and comfort of your EVOKE™, walk or perform several flexions.


Adjust the tension of the straps if necessary.



\* Some of these steps may not apply to your EVOKE™.



# EVOKE™



# OssKin

Revolution in movement


## A MESSAGE FROM THE TEAM

Welcome to the Osskin family!

The Osskin team has set itself the goal of making the wearing experience of your knee brace the most pleasant and comfortable possible, while also significantly reducing your pain.

After more than 15 years of R&D, numerous patents, and many clinical trials, **EVOKE™** was born.

That is why we are proud to offer you **EVOKE™**, the lightest and most discreet knee brace in the world.



## UNIQUE

Did you know that your **EVOKE™** was produced using a 3D printing manufacturing process that was developed in the aerospace industry?

All of this to make a highly custom product that perfectly respects your anatomical profile.

Working with this advanced technology allows us to obtain a maximum of geometric complexity, while ensuring a superior level of support and comfort and unmatched precision.

Your leg and your pain are unique. We want to offer you a product that is as well.

## CAUTION

- ◆ If you feel abnormal pain, change of sensations or unusual reactions while wearing the **EVOKE™**, contact your professional orthotist or consult your doctor.

## HINGE SYSTEM

- ◆ For safety reasons, the hinges must be disassembled only by your professional orthotist;
- ◆ Regularly check hinges for foreign objects after each use and remove anything found in the hinges of your **EVOKE™**;
- ◆ Clean hinges with clean water or compressed air;
- ◆ Never lubricate the hinges. If the hinges do not work normally, they must be repaired by your professional orthotist.

## CARE AND CLEANING


- ◆ Salt water, sand and all types of dirt can damage your product. If exposed to these conditions, rinse thoroughly with clean water.

## NOTE:

Maintain your **EVOKE™** in order to preserve the optimal condition of the frame, hinges, straps, and cushioning.

## TIPS AND ADVICE

- ◆ Apply your **EVOKE™** directly to the skin, leggings or tights. Do not wear your product over rigid pants.
- ◆ During sport or physical activities, wearing a legging or tights is recommended and must be worn under your **EVOKE™**.
- ◆ Do not apply too much tension to the femoral straps #2 and #3.



## WARRANTY

|               |          |
|---------------|----------|
| Frame/hinges  | 5 years  |
| Pads & liners | 6 months |
| Straps        | 6 months |

Osskin guarantees its products from the date of purchase against any material and manufacturing defects. Normal wear and tear during use is not considered a defect. If **EVOKE™** is modified by someone else other than your professional orthotist or your doctor, the warranty will no longer apply. Contact your professional orthotist for more details.

