










INSTRUCTIONS FOR APPLYING THE PROVOKE BRACE


- A**  **B** 
- **GETTING READY**
- Sit on the edge of a stable chair and bend your knee at a 10° angle (A). Attach the straps to themselves so they don't get in the way (Fig. B).

- C**  **D** 
- **ALIGNING THE BRACE**
- Align the center of the hinges (C) along the center of the condyles of your knee (D).

- E** 
- **STRAP #1**
- Attach the proximal posterior tibial strap #1 over the calf to keep the brace in place (E).

- F**  **G** 
- **STRAP #2 AND #3**
- Attach the posterior proximal femoral strap #2 (F) and then the posterior distal tibia strap #3 if present (G).

- H**  **I** 
- **STRAP #4 AND #5**
- Attach the posterior distal femoral strap #4 (H) and the anterior tibia strap #5, if present (I).

- J** 
- **ONE LAST ADJUSTMENT**
- Stand up and tighten the strap around your thigh (J). Check the fit and comfort of your PROVOKE, by walking or performing a few knee bends. Adjust the straps as necessary.

Provoke™



WELCOME TO OUR TEAM!

OSSKIN IS HAPPY AND PROUD THAT YOU'VE JOINED OUR FAMILY!

WHEN WE STARTED, WE SET OURSELVES THE GOAL OF CREATING THE MOST PLEASANT AND COMFORTABLE KNEE-BRACE EXPERIENCE POSSIBLE, WHILE MARKEDLY REDUCING YOUR PAIN AND HELPING YOU GET BACK TO LIVING!

FOLLOWING 15 YEARS OF INTENSIVE R&D, MANY PATENTS AND CLINICAL TRIALS, WE WERE PROUD AND HAPPY TO ANNOUNCE THE LAUNCH OF PROVOKE, THE WORLD'S LIGHTEST AND MOST DISCREET KNEE BRACE.

TRULY ONE OF A KIND

YOUR PROVOKE BRACE WAS CREATED USING A 3D PRINTING PROCESS DEVELOPED BY THE AEROSPACE INDUSTRY TO ENSURE A FULLY CUSTOMIZED PRODUCT THAT MATCHES YOUR LEG'S ANATOMY PERFECTLY.

USING THIS STATE-OF-THE-ART TECHNOLOGY ENSURES US MAXIMUM GEOMETRIC COMPLEXITY AND UNBEATABLE PRECISION FOR SUPERIOR SUPPORT AND COMFORT.

YOUR LEG AND THE PAIN YOU EXPERIENCE ARE UNIQUE TO YOU. WITH PROVOKE WE OFFER YOU A BRACE THAT'S EVERY BIT AS UNIQUE, ONE THAT WILL FOLLOW YOU EVERY STEP OF THE WAY!

ASYMMOTION HINGE SYSTEM

To ensure safety, the hinges must only be disassembled by a professional orthotist;

Check hinges after every use and remove any foreign objects or debris.

Clean hinges with clean water or compressed air;

Never lubricate the hinges.

If the hinges do not work normally, they must be repaired by your professional orthotist.

CLEANING AND CARE

Dirt, sand and salt water can all damage your brace. If your brace comes into contact with any of the above, rinse it thoroughly with clean water.

Hand wash with detergent-free soap.

Rinse well with cold water and allow to air dry.

NOTE: Straps, pads and liners are all detachable.

HELPFUL TIPS

Wear your PROVOKE directly on your skin, leggings or tights. Do not wear over rigid pants.

Wearing your PROVOKE over leggings or tights is highly recommended during sports or physical activities.

Do not apply too much tension to the femoral straps #2 and #3.

WARRANTY

› **FRAME/HINGES: 3 YEARS**

› **PADS & LINERS: 6 MONTHS**

› **STRAPS: 6 MONTHS**

OssKin guarantees its products from the date of purchase against any material and manufacturing defects. Normal wear and tear during use is not considered a defect. If PROVOKE is modified by someone else other than a professional orthotist or your doctor, the warranty will no longer apply. Contact your professional orthotist for more details.

CAUTION!

If you feel abnormal pain, change of sensations or unusual reactions while wearing the PROVOKE, contact your professional orthotist or consult your doctor.

